According to the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number. The valid OMB control number for this information collection is 0579-0315. The time required to complete this information collection is estimated to average 2 hours per response, including the time to review instructions, search existing data resources, gather the data needed, and complete and review the information collected.

OMB Approved 0579-0315 EXP: XX/20XX

UNITED STATES DEPARTMENT OF AGRICULTURE
ANIMAL AND PLANT HEALTH INSPECTION SERVICE
VETERINARY SERVICES
NATIONAL ANIMAL HEALTH MONITORING SYSTEM
2150 CENTRE AVE, BLDG B
FORT COLLINS, CO 80526

Swine 2020 Large Enterprise Saliva Collection Record

Kit Contents

- 30 feet of 5/8" 3 strand cotton rope
- 10 quart Ziploc® bags
- 1 gallon Ziploc® bag
- 10, 50 ml centrifuge tubes
- Insulated cooler and shipping box

- 3 ice packs (freeze prior to farm visit)
- 1 absorbent sheet
- 2 linear bags (one used as trash bag, one used to line the cooler)
- Paperwork packet: Oral Fluids Collection Record, 10 centrifuge tubes labels, shipping airbill

Collection Instructions

Collect samples on a **Sunday-Thursday**. It is important that **the site house 20 weeks and older finishing pigs in a barn or enclosure (i.e., a hoop barn) in pens**. Hang ropes in a maximum of **8 pens in one barn or room within a barn** in pens populated by these pigs. If one barn or room contains the site's oldest finishing pigs, pick that one. Otherwise take any barn that is available with the right age and type of pigs. The ropes should be hung in pens that are roughly equidistant from one another.

Ropes are hung by taking one end of the cotton rope (5/8" 3 Strand Cotton) and tying one end to the top bar of each pen or the gate. Tie the rope knot on the gate and cut the end in the pen so that that this end hangs approximately at shoulder height for most of the animals in the pen. Unravel the part of the rope that hangs in the pen rope into its strands. Shake the end of the rope in the pen so that pigs are attracted to it.

Wait 30-45 minutes while they chew on the ropes. While the rope is still hanging, insert the lower, wet portion of each rope into a quart sized Ziploc® bag and cut the part that's inside the bag from the remainder of the hanging rope. Manually squeeze the rope pieces while inside the Ziploc® bag to release the oral fluid. Cut the corner of the Ziploc® bag and let the fluid drip into a 50 ml centrifuge tube. Close each centrifuge tube, checking that the lid is tight. Place the used Ziploc® bags and the rope ends in the trash bag. Untie and place excess hanging rope in the trash bag. Write the **8-digit Farm ID** (State/operation/site), **Facility, and Pen IDs on a sample label,** peel off and place securely on the appropriate centrifuge tubes.

Place centrifuge tubes in the gallon sized Ziploc® plastic bag and seal. Place the Ziploc® bag inside the liner bag with 3 frozen ice packs and the absorbent sheet. Place the liner bag inside the insulated shipment container and close the shipping container lid.

Collection Record Instructions

Complete the Oral Fluids Collection Record Form (attached below). Place the original (white) collection record on top of the insulated cooler box lid. Leave the yellow copy of the collection record with the producer. Secure the shipping box.

Shipping and Results

Keep samples cool in a refrigerator, attach shipping airbill (addressed to ISU) to the shipping box, and ship via **Fed Ex within 24 hours of collection on Monday-Thursday**. The samples will be tested for **PRRS**, **Flu**, **SVA and three types of Coronaviruses** and results will be returned to the Producer in approximately 45-60 days after collection.

Swine 2020: Oral Fluids Collection Record

Farm ID	# of People Involved	Key Collector	Collection Date	Hours to Take	Hours of Overall	Kit Number
	Fed VMO Fed AHT	Initials		and Prep	Travel Time	
	Producer State VMO			Samples		
	St AHT Priv Vet					
	Others-specify:					
8-digits: St, Op, Site	Enter number for each category		MM/DD/YY	In quarter hours	In quarter hours	Must match labels

Rope			Facility	Where was	Approximate # of pigs in	# of pigs share air	Average age of pigs in pen		Vac	Vaccination codes ⁶		
	Facility/Barn ID ¹	Pen ID ¹	•	rope hung?3		space ⁴	(weeks)	Gender of pigs in pen5	1st	2nd	3rd	4th
1			□T □O	□B □G □O				□G □B □M				
2			□T □O	□B □G □O				□G □B □M				
3			□T □O	□B □G □O				□G □B □M				
4			□T □O	□B □G □O				□G □B □M				
5			□T □O	□B □G □O				□G □B □M				
6			□T □O	□B □G □O				□G □B □M				
7			□T □O	□B □G □O				□G □B □M				
8			□T □O	□B □G □O				□G □B □M				

¹Any designation will do whether text or numeric as long as it differentiates between row entries.

⁶Up to **4** different vaccinations may be entered. The codes and what they stand for are listed below. These are vaccinations that may have been given at any time in the pigs' lives.

Vaccination codes				
1=APP (Actinobacillus pleuropneumoniae)	11=Influenza			
2=Actinobacillus suis (autogenous)	12=Leptospirosis			
3=Atrophic rhinitis (Bordatella/Pasteurella)	13=Mycoplasma hyopneumoniae			
4=Clostridium difficile (autogenous)	14=Porcine circovirus 2			
5=Clostridium perfringens Type A	15=PRRS			
6=Clostridium perfringens Types C and D	16=Porcine epidemic diarrhea			
7=Erysipelas	17=Rotavirus			
8=E. coli (K88, K99, 987P, F41)	18=Salmonella			
9=Glasser's disease (Haemophilus parasuis)	19=Streptococcus suis			
10=lleitis (Lawsonia intracellularis)	20=TGE (transmissible gastroenteritis)			

²T=Total Confinement, O=Open building. Total confinement means solid exterior walls and mechanical ventilation. An open building generally has open sides or the potential for sides to be opened weather permitting.

³B=Bar (e.g., top bar of pen enclosure), G=Gate, O=Other

⁴This would be the number of pigs in the barn or room within a barn usually.

⁵G=Gilts, B=Barrows and M=Mixed.