

At-Home OTC COVID-19 Diagnostic Tests

Español (/medical-devices/coronavirus-covid-19-and-medical-devices/pruebas-de-covid-19-sin-receta-medica-o-de-venta-libre-para-diagnostico-en-el-hogar)



On this page:

- [General Information About At-Home OTC COVID-19 Diagnostic Tests](#)
- [List of Authorized At-Home OTC COVID-19 Diagnostic Tests](#)
- [Related Information](#)

General Information About At-Home OTC COVID-19 Diagnostic Tests

This page has information about at-home over-the-counter (OTC) COVID-19 diagnostic tests. Diagnostic tests can show if you have an active COVID-19 infection.

These at-home OTC COVID-19 diagnostic tests are [FDA authorized \(https://www.fda.gov/emergency-preparedness-and-response/mcm-legal-regulatory-and-policy-framework/emergency-use-authorization\)](https://www.fda.gov/emergency-preparedness-and-response/mcm-legal-regulatory-and-policy-framework/emergency-use-authorization) for self-testing at home (or in other locations). This means you collect your own sample, perform the test, and read the result yourself without the need to send a sample to a laboratory.

There are two types of COVID-19 diagnostic tests, molecular and antigen, which detect different parts of the virus. Most authorized at-home OTC COVID-19 tests are antigen tests, and there are also a small number of authorized at-home OTC COVID-19 molecular tests.

When using a COVID-19 antigen diagnostic test, you should repeat testing following a negative result, whether you have symptoms or not, to reduce your risk of a false negative test result. For more information about this, read our [FDA Safety Communication \(https://www.fda.gov/medical-devices/safety-communications/home-covid-19-antigen-tests-take-steps-reduce-your-risk-false-negative-fda-safety-communication\)](https://www.fda.gov/medical-devices/safety-communications/home-covid-19-antigen-tests-take-steps-reduce-your-risk-false-negative-fda-safety-communication).

For additional information on reading and understanding your test results, see [Understanding At-Home OTC COVID-19 Antigen Diagnostic Test Results \(/medical-devices/coronavirus-covid-19-and-medical-devices/understanding-home-otc-covid-19-antigen-diagnostic-test-results\)](/medical-devices/coronavirus-covid-19-and-medical-devices/understanding-home-otc-covid-19-antigen-diagnostic-test-results).

Authorized at-home OTC tests are available without a prescription and are typically available online or at local stores. Testing and other protective steps like wearing a mask and COVID-19 vaccination are important to stop the spread of COVID-19 infection.

For more information, including answers to frequently asked general questions about at-home COVID-19 diagnostic tests, visit [At-Home COVID-19 Diagnostic Tests: Frequently Asked Questions \(/medical-devices/coronavirus-covid-19-and-medical-devices/home-covid-19-diagnostic-tests-frequently-asked-questions\)](/medical-devices/coronavirus-covid-19-and-medical-devices/home-covid-19-diagnostic-tests-frequently-asked-questions).

List of Authorized At-Home OTC COVID-19 Diagnostic Tests

The table below lists FDA-authorized at-home OTC COVID-19 diagnostic tests, and includes information on expiration dates, who can use the test, and other details that may help you decide what test is right for you. To see complete information on smaller screens, select the blue plus (+) button beside the test name.

The table includes links to home use instructions for each test. For more information about each test, including the Letter of Authorization and authorized labeling, see [In Vitro Diagnostics EUAs: Tables of IVD EUAs \(/medical-devices/coronavirus-disease-2019-covid-19-emergency-use-authorizations-medical-devices/in-vitro-diagnostics-euas\)](/medical-devices/coronavirus-disease-2019-covid-19-emergency-use-authorizations-medical-devices/in-vitro-diagnostics-euas).

In the table below, the "Expiration Date" column lists the shelf-life and where to find the expiration date for that test. The shelf-life is how long the test should work as expected and is measured from the date the test was manufactured. The expiration date is set at the end of the shelf-life. In some cases, the expiration date for a test may be extended. An extended expiration date means the manufacturer provided data showing that the shelf-life is longer than was known when the test was first authorized. For more information about how the expiration date is determined and why it may be extended, see the [At-Home COVID-19 Diagnostic Tests: Frequently Asked Questions \(/medical-devices/coronavirus-covid-19-and-medical-devices/home-covid-19-diagnostic-tests-frequently-asked-questions\)](/medical-devices/coronavirus-covid-19-and-medical-devices/home-covid-19-diagnostic-tests-frequently-asked-questions).

To see if the expiration date for your at-home OTC COVID-19 test has been extended, first find the row in the below table that matches the manufacturer and test name shown on the box label of your test.



- If the Expiration Date column says that the shelf-life is "extended," there is a link to "updated expiration dates" where you can find a list of the original expiration dates and the new expiration dates. Find the original expiration date on the box label of your test and then look for the new expiration date in the "updated expiration dates" table for your test.
- If the Expiration Date column **does not** say the shelf-life is extended, that means the expiration date on the box label of your test is still correct. The table will say "See box label" instead of having a link to updated expiration dates.





The FDA will update this table as additional shelf-life extensions are authorized.

To see complete information on smaller screens, select the blue plus (+) button beside the test name.

Search:



Show entries



Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
<p>Abbott Diagnostics Scarborough, Inc.: BinaxNOW COVID-19 Ag Card Home Test (/media/144575/download)</p>  <p>(/media/160087/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 15 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Requires supervision of a telehealth proctor and a smartphone or computer • Results in 15 minutes 	<ul style="list-style-type: none"> • Shelf Life: 15 months (extended) • Expiration Date: See link for list of updated expiration dates (/media/158002/download)
<p>Abbott Diagnostics Scarborough, Inc.: BinaxNOW COVID-19 Antigen Self Test (/media/147255/download)</p>  <p>(/media/160094/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 15 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Smartphone optional • Results in 15 minutes 	<ul style="list-style-type: none"> • Shelf Life: 22 months (extended) • Expiration Date: See link for list of updated expiration dates (/media/158003/download)

Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
<p>Access Bio, Inc.: CareStart COVID-19 Antigen Home Test (/media/152188/download)</p>  <p>(/media/160095/download)</p> <p>Alternate brand name:</p>  <p>(/media/160097/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Smartphone optional • Alternate brand name: On/go COVID-19 Antigen Self-Test • Results in 10 minutes 	<ul style="list-style-type: none"> • Shelf Life: 15 months (extended) • Expiration Date: See link for list of updated expiration dates (/media/158734/download)
<p>ACON Laboratories, Inc: Flowflex COVID-19 Antigen Home Test (/media/152699/download)</p>  <p>(/media/160099/download)</p> <p>Alternate brand name:</p>  <p>(/media/160098/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Smartphone optional • Alternate brand name: On/Go One COVID-19 Antigen Home Test • Results in 15 minutes 	<ul style="list-style-type: none"> • Shelf Life: 21 months (extended) • Expiration Date: See link for list of updated expiration dates (/media/158005/download)


Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
Advin Biotech Inc.: Advin COVID-19 Antigen Test @Home (/media/164171/download)	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days. The test is to be performed two times over three days (serial testing). • People without symptoms. The test is to be performed three times over five days (serial testing). 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 15-30 minutes 	<ul style="list-style-type: none"> • Shelf Life: 12 months • Expiration Date: See box label






Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
<p>ANP Technologies, Inc.: NIDS COVID-19 Antigen Home Test (/media/163321/download).</p>  <p>(/media/163324/download).</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days. The test is to be performed two times over three days (serial testing). • People without symptoms. The test is to be performed three times over five days (serial testing). 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 15-30 minutes 	<ul style="list-style-type: none"> • Shelf Life: 11 months • Expiration Date: See box label
<p>Aptitude Medical Systems Inc.: Metrix COVID-19 Test (https://www.fda.gov/media/162404/download).</p>  <p>(/media/162438/download).</p>	<ul style="list-style-type: none"> • People with or without symptoms 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Molecular test • Nasal swab • Saliva • Requires a Metrix Reader (sold separately) • Smartphone optional • Results in 30 minutes 	<ul style="list-style-type: none"> • Shelf Life: 12 Months • Expiration Date: See box label




Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
<p>Azure Biotech Inc.: Fastep COVID-19 Antigen Home Test (/media/163448/download)</p>  <p>(/media/163451/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 5 days. The test is to be performed two times over three days (serial testing). • People without symptoms. The test is to be performed three times over five days (serial testing). 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 15-30 minutes 	<ul style="list-style-type: none"> • Shelf Life: 24 months • Expiration Date: See box label
<p>Becton, Dickinson and Company (BD): BD Veritor At-Home COVID-19 Test (/media/151763/download)</p>  <p>(/media/160100/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Requires a smartphone • Results in 15 minutes 	<ul style="list-style-type: none"> • Shelf Life: 6 months • Expiration Date: See box label

Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
<p>Beijing Hotgen Biotech Co., Ltd.: Hotgen COVID-19 Antigen Home Test (/media/163439/download)</p>  <p>(/media/163364/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days. The test is to be performed two times over three days (serial testing). • People without symptoms. The test is to be performed three times over five days (serial testing). 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 15-30 minutes 	<ul style="list-style-type: none"> • Shelf Life: 6 months • Expiration Date: See box label
<p>Celltrion USA, Inc.: Celltrion DiaTrust COVID-19 Ag Home Test (/media/153421/download)</p>  <p>(/media/160101/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older 	<ul style="list-style-type: none"> • Antigen test • Mid-turbinate nasal swab • Smartphone optional • Results in 15 minutes 	<ul style="list-style-type: none"> • Shelf Life: 18 months (extended) • Expiration Date: See link for updated expiration dates (/media/159236/download)



Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
<p>CorDX, Inc.: CorDx COVID-19 Ag Test (https://www.fda.gov/media/164118/download)</p>  <p>(/media/163450/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days. The test is to be performed two times over three days (serial testing). • People without symptoms. The test is to be performed three times over five days (serial testing). 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 10-30 minutes 	<ul style="list-style-type: none"> • Shelf Life: 15 months • Expiration Date: See box label




Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
CTK Biotech, Inc.: ImmuView COVID-19 Antigen Home Test (/media/164140/download) 	<ul style="list-style-type: none"> • People with symptoms that began within the last 5 days. The test is to be performed two times over three days (serial testing). • People without symptoms. The test is to be performed three times over five days (serial testing). 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 15-30 minutes 	<ul style="list-style-type: none"> • Shelf Life: 12 months • Expiration Date: See box label
Cue Health Inc.: Cue COVID-19 Test for Home and Over The Counter (OTC) Use (/media/146471/download)  (/media/160102/download)	<ul style="list-style-type: none"> • People with or without symptoms 	<ul style="list-style-type: none"> • Age 18 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Molecular test • Nasal swab • Requires a Cue Cartridge Reader (sold separately) • Requires a smartphone • Results in 20 minutes 	<ul style="list-style-type: none"> • Shelf Life: 4 months • Expiration Date: See box label





Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
<p>Detect, Inc.: Detect Covid-19 Test (/media/153746/download)</p>  <p>(/media/160103/download)</p>	<ul style="list-style-type: none"> • People suspected of COVID-19 • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Molecular test • Nasal swab • Requires a Detect Hub (sold separately) • Smartphone optional • Results in 60 minutes 	<ul style="list-style-type: none"> • Shelf-Life: 8 months (extended) • Expiration Date: See link for list of updated expiration dates (/media/158004/download)
<p>Ellume Limited: Ellume COVID-19 Home Test (/media/144593/download)</p>  <p>(/media/160104/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 16 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Mid-turbinate nasal swab • Requires a Smartphone • Results in 15 minutes 	<ul style="list-style-type: none"> • Shelf Life: 12 months • Expiration Date: See box label



Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
<p>Genabio Diagnostics Inc.: Genabio COVID-19 Rapid Self-Test Kit (/media/159834/download)</p>  <p>(/media/160105/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days. • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 15-30 minutes 	<ul style="list-style-type: none"> • Shelf Life: 18 months • Expiration Date: See box label
<p>iHealth Labs, Inc.: iHealth COVID-19 Antigen Rapid Test (/media/153924/download)</p>  <p>(/media/160107/download)</p> <p>Alternate brand name:</p>  <p>(/media/160106/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 15 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Smartphone optional • Results in 15 minutes • Alternate brand name: GoodToKnow COVID-19 Antigen Rapid Test • Alternate brand name: GoToKnow COVID-19 Antigen Rapid Test 	<ul style="list-style-type: none"> • Shelf Life: 15 months (extended) • Expiration Date: See link for list of updated expiration dates (/media/164551/download)


Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
<p>InBios International Inc: SCoV-2 Ag Detect Rapid Self-Test (/media/154337/download)</p>  <p>(/media/160108/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 5 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 20 minutes 	<ul style="list-style-type: none"> • Shelf Life: 13 months (extended) • Expiration Date: See link for list of updated expiration dates (/media/158008/download)
<p>Lucira Health, Inc: Lucira CHECK-IT COVID-19 Test Kit (/media/147495/download)</p>  <p>(/media/160109/download)</p>	<ul style="list-style-type: none"> • People with or without symptoms 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Molecular test • Nasal swab • Smartphone optional • Results in 30 minutes 	<ul style="list-style-type: none"> • Shelf Life: 18 months (extended) • Expiration Date: See link for list of updated expiration dates (/media/161056/download)
<p>Maxim Biomedical, Inc.: MaximBio ClearDetect COVID-19 Antigen Home Test (/media/155635/download)</p>  <p>(/media/160110/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 5 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 15 minutes 	<ul style="list-style-type: none"> • Shelf Life: 12 months (extended) • Expiration Date: See link for list of updated expiration dates (/media/160581/download)

Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
<p>Oceanit Foundry LLC: ASSURE-100 Rapid COVID-19 (/media/164167/download)</p> 	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days. The test is to be performed two times over three days (serial testing). • People without symptoms. The test is to be performed three times over five days (serial testing). 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 20-30 minutes 	<ul style="list-style-type: none"> • Shelf Life: 12 months • Expiration Date: See box label
<p>OraSure Technologies, Inc.: InteliSwab COVID-19 Rapid Test (/media/149912/download)</p>  <p>(/media/160111/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 18 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Smartphone optional • Results in 30 minutes 	<ul style="list-style-type: none"> • Shelf Life: 12 months (extended) • Expiration Date: See link for list of updated expiration dates (/media/161057/download)

Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
<p>OSANG LLC: OHC COVID-19 Antigen Self Test (/media/157549/download)</p>  <p>(/media/161361/download) Alternate brand name:  (/media/161362/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days. The test is to be performed two times over three days (serial testing). • People without symptoms. The test is to be performed two times over three days (serial testing). 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Alternate brand name: QuickFinder COVID-19 Antigen Self Test • Results in 15 minutes 	<ul style="list-style-type: none"> • Shelf Life: 12 months (extended) • Expiration Date: See link for list of updated expiration dates (/media/161358/download)
<p>PHASE Scientific International, Ltd.: INDICAID COVID-19 Rapid Antigen At-Home Test (/media/156955/download)</p>  <p>(/media/160113/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 6 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 20 minutes 	<ul style="list-style-type: none"> • Shelf Life: 12 months • Expiration Date: See box label

Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
<p>Quidel Corporation: QuickVue At-Home OTC COVID-19 Test (/media/147250/download)</p>  <p>(/media/162058/download)</p> <p>Alternative Brand Names:</p>  	<ul style="list-style-type: none"> • People with symptoms that began within the last 6 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Alternative Brand Names: CVS Health At Home COVID-19 Test Kit and Walgreens At-Home COVID-19 Test Kit • Results in 10 minutes 	<ul style="list-style-type: none"> • Shelf Life: 16 months (extended) • Expiration Date: See link for updated expiration date (/media/162056/download)
<p>SD Biosensor, Inc.: Pilot COVID-19 At-Home Test (/media/155126/download)¹</p>  <p>(/media/160115/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 6 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 20 minutes • Telehealth Proctor Supervised (optional) 	<ul style="list-style-type: none"> • Shelf Life: 12 months (extended) • Expiration Date: See link for updated expiration date (/media/159640/download)

Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
<p>Siemens Healthineers: CLINITEST Rapid COVID-19 Antigen Self-Test (/media/155176/download)</p>  <p>(/media/160116/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 15 minutes 	<ul style="list-style-type: none"> • Shelf Life: 24 months (extended) • Expiration Date: See link for updated expiration date (/media/162503/download)
<p>Watmind USA: Speedy Swab Rapid COVID-19 Antigen Self-Test (/media/159828/download)</p>  <p>(/media/160117/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 6 days. The test is to be performed two times over three days (serial testing) • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 15 minutes 	<ul style="list-style-type: none"> • Shelf Life: 6 months • Expiration Date: See box label

Manufacturer and Test Name (Links to Instructions for Use)	Who can use this test: Symptoms	Who can use this test: Age	Other Details	Expiration Date
<p>Xiamen Boson Biotech Co., Ltd.: Rapid SARS-CoV-2 Antigen Test Card (/media/157545/download)</p>  <p>(/media/160118/download)</p>	<ul style="list-style-type: none"> • People with symptoms that began within the last 6 days. The test is to be performed two times over three days (serial testing). • People without symptoms. The test is to be performed two times over three days (serial testing). 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Antigen test • Nasal swab • Results in 15 minutes 	<ul style="list-style-type: none"> • Shelf Life: 6 months • Expiration Date: See box label

Showing 1 to 30 of 30 entries

Previous **1** Next

Related Information:

- [Understanding At-Home OTC COVID-19 Antigen Diagnostic Test Results](#) ([/medical-devices/coronavirus-covid-19-and-medical-devices/understanding-home-otc-covid-19-antigen-diagnostic-test-results](#))
- [At-Home COVID-19 Diagnostic Tests: Frequently Asked Questions](#) ([/medical-devices/coronavirus-covid-19-and-medical-devices/home-covid-19-diagnostic-tests-frequently-asked-questions](#))
- [COVID-19 Test Basics:](#) ([/consumers/consumer-updates/covid-19-test-basics](#)) Includes details on COVID-19 tests, types of samples, and other information.
- [At-Home COVID-19 Antigen Tests-Take Steps to Reduce Your Risk of False Negative: FDA Safety Communication](#) ([/medical-devices/safety-communications/home-covid-19-antigen-tests-take-steps-reduce-your-risk-false-negative-results-fda-safety](#))
- [Counterfeit At-Home OTC COVID-19 Diagnostic Tests](#) (<https://www.fda.gov/medical-devices/coronavirus-covid-19-and-medical-devices/counterfeit-home-otc-covid-19-diagnostic-tests>)
- Find COVID-19 Tests at [COVID.gov/tests](https://www.covid.gov/tests) (<https://www.covid.gov/tests>)
- Video: [How to Interpret Self-Test Results | CDC](https://www.youtube.com/watch?v=yoZL1BVP9Bs) (<https://www.youtube.com/watch?v=yoZL1BVP9Bs>) [↗](#)
(<http://www.fda.gov/about-fda/website-policies/website-disclaimer>)

- [COVID-19 Testing: What You Need to Know | CDC \(https://www.cdc.gov/coronavirus/2019-ncov/symptoms-testing/testing.html\)](https://www.cdc.gov/coronavirus/2019-ncov/symptoms-testing/testing.html)
- [COVID-19 Self-Testing At-Home or Anywhere | CDC \(https://www.cdc.gov/coronavirus/2019-ncov/testing/self-testing.html\)](https://www.cdc.gov/coronavirus/2019-ncov/testing/self-testing.html)
- [What to Do If You Were Exposed to COVID-19 | CDC \(https://www.cdc.gov/coronavirus/2019-ncov/your-health/if-you-were-exposed.html\)](https://www.cdc.gov/coronavirus/2019-ncov/your-health/if-you-were-exposed.html)
- [Isolation and Precautions for People with COVID-19 | CDC \(https://www.cdc.gov/coronavirus/2019-ncov/your-health/isolation.html\)](https://www.cdc.gov/coronavirus/2019-ncov/your-health/isolation.html)
- Video: [What Is an EUA? \(https://www.youtube.com/watch?app=desktop&v=iGkwaESsGBQ&feature=youtu.be\)](https://www.youtube.com/watch?app=desktop&v=iGkwaESsGBQ&feature=youtu.be)
[↗ \(http://www.fda.gov/about-fda/website-policies/website-disclaimer\)](http://www.fda.gov/about-fda/website-policies/website-disclaimer) Describes how the FDA can issue an emergency use authorization (EUA) to provide more timely access to diagnostic tests that may help during the public health emergency when there is no adequate, approved, and available alternative.

¹ The Pilot COVID-19 At-Home Test, distributed by Roche Diagnostics, is previously known as the COVID-19 At-Home Test. As the manufacturer, SD Biosensor, transitions to this new brand, customers may receive a product with the previously authorized COVID-19 At-Home Test labeling.