



Vitamin A Laboratory - External Quality Assurance Program (VITAL-EQA)
Round 35 / Fat Solubles / Retinol

Form Approved
 OMB No. 0920-xxxx
 Exp. Date xx/xx/20xx

INSTRUCTIONS

- E-mail this completed worksheet to: **vitaminalab@cdc.gov**
- Round deadline: **December 18, 2020**
- All cells highlighted **red** are required or have invalid data that needs to be revised. Click on cell to see data requirements.
- Run these samples in the same way as you would routine patient samples.
- Run in singlicate over a period of two consecutive days (Day 1 and Day 2).
- Record all results with three significant figures (i.e., 105, 10.5, 1.05, 0.105).
- Retain the vials in your ultra-cold freezer for at least two weeks after data submission in case there is a question about the ID code or results.

Laboratory ID

Received Date (MM/DD/YY) res, snipper arrived cold no, snipper did not arrive cold

Assay Type

Instrument (Manufacturer/Model)
 Kit Information (if used) (Name/Manufacturer)

Calibration Range Low High $\mu\text{g/dL}$ $\mu\text{g/dL}$

LOD $\mu\text{g/dL}$

Run	Assay Date	Sample ID	Result ($\mu\text{g/dL}$)	< LOD	Note
1			1	<input type="checkbox"/>	
			1	<input type="checkbox"/>	
			1	<input type="checkbox"/>	
2			1	<input type="checkbox"/>	
			1	<input type="checkbox"/>	
			1	<input type="checkbox"/>	

Laboratory Notes

OFFICIAL USE ONLY